

Care Vision DA10

Automated External Defibrillator



Features

- Support Icons and voice prompts to instruct user performing CPR
- By observing the visual indicator LEDs, easy following of operating steps
- An automatic self-test which is performed on a daily, weekly and monthly
- After analyzing ECG, flashing Shock button LED indicates clearly readiness for shock delivery
- Using SD card, customer can check ECG status from PC or upgrade the unit



Care Vision Health Care Pte. Ltd.

DA10

Care Vision DA10

Automated External Defibrillator

Specification

Dimensions	314x259x109 (mm, HxWxD)
Weight	Approx. 2.95kg(with battery)
Defibrillation	Waveform :BiphasicSmartshot, Energy: 45~200J, Charging Time: max 13 sec, Operating mode: Semi-Auto
ECG	lead: II, VF/VT detection, Analysis time: less than 13 sec.
Indication	Icon, Status LCD, Voice prompt, CPR indication
Data backup &	SD card, IrDA, LAN(option)
Battery	Type: LiMnO ₂ , Storage: 4 years, Discharge: Min 200, Monitoring: Min 10hrs
Regulatory	AAMI DF80, IEC 60601-1-2, IEC 60601-1-4, IEC 60601-2-4, 2010 AHA guideline
Self test	Every 24hr, 1wek, 1month, Power on self test, Test result: Status LCD
PADS	Adult pad: 1 disposable pad (standard Acc.) - Storage: 2 years or Pediatric Pad: 1 disposable pad (standard Acc.) - Storage: 2 years
List of Components	DA10, Adult pad, Operation Manual



IrDA port

IR communication port between PC

SD Card

Review date stored & software upgrade



Battery

Disposable LiMnO₂ (non-rechargeable)



Adult Pads



Pediatric Pads



* Technical specifications are subject to change without prior notice. The functions and technical features detailed here are valid at the time of press. They are neither legally binding nor guaranteed. All features are liable to change without notice. For latest data, please contact your dealer.